



A step ahead

FujiCEM™ Evolve  
from GC

Innovative resin-modified glass  
ionomer cement for  
a trustful performance and  
enjoyable convenience

**GC**

# FujiCEM Evolve from GC

## A step ahead

Glass ionomers are GC's legacy, with proven quality as seen in independent studies all over the world. Take a step ahead by choosing for easy application, less technique sensitivity and high moisture tolerance with FujiCEM Evolve.



### Hand mix or automix



- Openings are far apart: pastes do not contaminate each other, hand mixing is possible

- Easy to place  
"Push and Click" mixing tip
- Can be dispensed directly into the restoration or the root canal (with the endo tip)



### The resin-modified glass ionomer you can trust

- Less technique-sensitive than resin cements
- No need for pre-treatment: one-step application
- Long working time, short setting time
- Moisture tolerant: rubber dam is not required
- Less risk of post-operative sensitivity
- Excellent sealing and marginal adaptation



Courtesy of Dr. J. Tapia Guadix, Spain

### Optimal viscosity


- Controlled extrusion for an easier procedure
- Low film thickness & precise adaptation
- Tack-cure option for easier excess removal



Courtesy of Dr. J. Tapia Guadix, Spain




## Reliable bonding in all circumstances



Source: GC R&D Data, Japan, 2018. Data on file.


## Optimised adhesion to zirconia



Source: GC R&D Data, Japan, 2018. Data on file.

## Long-chain monomer technology

- Super-long monomer chains for strength and stability of the cement.



- Radiopaque glass filler for excellent visibility on X-ray (> 250% AI)

## A versatile option for the everyday practice

FujiCEM Evolve is the versatile option for your everyday practice as it can be used for many indications and does not require absolute moisture isolation. Moreover, its improved bond strength to zirconia makes it particularly interesting to lute your high-strength all-ceramic restorations.

Indicated for:



# FujiCEM Evolve from GC

## A step ahead



Initial situation



Dispensing with automix



Seating of the restoration



Rubbery state is reached after just 3 s of tack-curing



Easy to remove excess at once



Final result

Courtesy of Prof. R. Sorrentino, Italy



### FujiCEM Evolve GC syringe (9.2 g / 5.0 ml)

012948	GC FujiCEM Evolve Starter Kit: Syringe (1), Tips Regular (15)
012949	GC FujiCEM Evolve Triple pack Handmix: Syringe (3), Mixing Pad (1)
012950	GC FujiCEM Evolve Triple pack Automix: Syringe (3), Tips Regular (45)
012954	GC Push and Click Tip Regular (15)
012955	GC Push and Click Tip Endo (15)



Find the step-by-step  
for all luting procedures  
on the **GC Luting Guide**



FujiCEM is a trademark of GC.

Ketac CEM Plus Luting, Meron Plus QM and Riva Cem are not trademarks of GC.

#### GC EUROPE N.V.

Head Office  
Researchpark  
Haasrode-Leuven 1240  
Interleuvenlaan 33  
B-3001 Leuven  
Tel. +32.16.74.10.00  
Fax. +32.16.40.48.32  
info.gce@gc.dental  
<http://www.gceurope.com>

#### GC UNITED KINGDOM Ltd.

Coopers Court  
Newport Pagnell  
UK-Bucks. MK16 8JS  
Tel. +44.1908.218.999  
Fax. +44.1908.218.900  
info.uk@gc.dental  
<http://uk.gceurope.com>

