



Simply The ONE

Reliable

Universal

Esthetic

Technique
insensitive

Simple



G-CEM ONE™
Self-adhesive resin cement



Since 1921
100 years of Quality in Dental

The ONE... that simplifies!

G-CEM ONE™ is a truly **universal, non-technique sensitive** cement with **high bond strength** and excellent self-curing ability. It simplifies **all your cementation procedures** and gives you more **predictable results** and peace of mind. G-CEM ONE™ can be used for a **wide range of indications**, including challenging clinical situations thanks to the **GC Touch Cure**.

GC Touch Cure

The chemical initiator in **G-CEM ONE™ ADHESIVE ENHANCING PRIMER** will accelerate the chemical cure of the cement from the tooth surface to ensure stable bonding.



↑ Touch Cure

G-CEM ONE™ ADHESIVE ENHANCING PRIMER

• Chemical Accelerator


Reliable

Optimal bond strength for retentive and non-retentive¹ preparations thanks to the unique **GC Touch Cure** ability of the optional G-CEM ONE™ ADHESIVE ENHANCING PRIMER.

Technique insensitive

Exceptional self curing for your peace of mind, tack cure feature for easy excess removal and clean up, moisture and saliva tolerance¹, optimized flow of the paste and virtually no post-operative sensitivity³.

One optimized self-curing ability



Source:
Prof. Miyazaki, Nihon University, Japan.

¹ With G-CEM ONE™ ADHESIVE ENHANCING PRIMER

² With G-Multi PRIMER™

³ Field survey results. Using G-CEM ONE™ only. Data on File.



The ONE... that simplifies!

G-CEM ONE™ ADHESIVE ENHANCING PRIMER (Optional)

For all abutment tooth and cavities

- MDP: bond to tooth structure, zirconia and non-precious metal
- 4-MET: bond to tooth structure
- MDTP: bond to precious metal
- Accelerator ⇒ Strong bonding to abutment tooth and cavity by **GC Touch Cure** feature.

No light curing, simply drying.

Wait 10 sec

Strong air 5 sec



Universal

Excellent bond strength to tooth and to all substrates:

Enamel, Dentin, Zirconia, Metal, Glass ceramics², Hybrid ceramics².

Esthetic

Invisible, wear resistant margin
for esthetic outcomes.

Simple

Reduced inventory for reproducibility
and optimized workflow.

One performance
regardless of the translucency



Metal




Lithium disilicate

Courtesy of
Dr. Kazunori Otani, Japan


Courtesy of
Dr. Akimasa Tsujimoto, Japan

The ONE... you can trust!

One cement for superior bond strength to all substrates




Source: GC R&D, Data on file.




Source: GC R&D, Data on file.
† With G-Multi PRIMER™

One primer for clinical challenges



Source: GC R&D, Data on file.


Immediate bond strength to avoid early debonding in non-retentive cases and CAD/CAM milled crowns, for the comfort of your patients.



Source: GC R&D, Data on file.

The G-CEM ONE™ ADHESIVE ENHANCING PRIMER ensures high bond strength, even in presence of saliva contamination.

One invisible margin for esthetic outcome



Source: GC R&D, Data on file.

Low film thickness (5 μm) and high wear resistance keep the cement line invisible.



Before treatment



After six months

Courtesy of Dr. Masayuki Okawa, Japan



With optional G-CEM ONE™ ADHESIVE ENHANCING PRIMER

The ONE... clinical partner!

One optimized flow



Ideal viscosity

One easy excess removal with 1 second tack curing



Courtesy of
Dr. Kazunori Otani, Japan

One for all Clinical case scenarios

Posterior Crown (Zirconia)



Prepared tooth cleaned and dried.



Restoration cemented.



Final result.

Courtesy of
Dr. Kazunori Otani, Japan



Anterior Crown (Lithium Disilicate)



Fractured veneer.



Crown prep.



Final result.

Dr. Miles Cone
Cerapist: Peter Pizzi, MDT



Posterior Onlay (Lithium Disilicate)



Initial situation.



Restoration cemented with use of
G-CEM ONE™ ADHESIVE ENHANCING
PRIMER.



Final result.

Courtesy of
Dr. Riccardo Tempesta, Italy



With optional G-CEM ONE™ ADHESIVE ENHANCING PRIMER

Packaging

SKU	Description
013670	G-CEM ONE™ Starter Kit A2
013671	G-CEM ONE™ Starter Kit Translucent
013672	G-CEM ONE™ System Kit
013673	G-CEM ONE™ Twin Refill A2
013674	G-CEM ONE™ Twin Refill Translucent
013675	G-CEM ONE™ Single Refill AO3
013676	G-CEM ONE™ Single Refill White Opaque
013642	G-CEM ONE™ ADHESIVE ENHANCING PRIMER



G-CEM ONE™ System Kit




G-CEM ONE™ Starter Kit



G-CEM ONE™ ADHESIVE ENHANCING PRIMER

Shade range



Related products



G-Multi PRIMER™
Universal Primer



GC Initial® LiSi Press
Lithium Disilicate



The Light
Curing Light



GC TEMP ADVANTAGE®
Temporary Cement

SKU #603125.04.2021.Ver1 • 800.323.7063 • www.gcamerica.com • www.gcamerica.com/training • © 2021 GC America Inc.

¹RelyX™ is a trademark of 3M, Maxcem Elite™ and NX3 Nexus™ are trademarks of Kerr Corporation, PANAVIA™ is a trademark of Kuraray Noritake Dental Inc., SpeedCem® and Multilink® are registered trademarks of Ivoclar Vivadent.

